

Transfer Turntable

The Turntable helps when transferring clients from one seated position to another.

The client stands up onto the Turntable which rotates smoothly through 360° to the new position.

Wipe clean finish.

Diameter: 16" (41 cm.)

Maximum User Weight: 19 ½ stone (125 Kg.)



Patient Hand Blocks

An aid to moving oneself in bed by raising the body above the mattress. Can be used singly or in pairs. Made from high strength polymer with a comfortable hand grip and non slip base pads.

Wipe clean finish. **CE Class I Medical Device.**

Height: 4" (10.5 cm.)

Base: 4 ¼ " x 8 ½ " (11 cm. x 22 cm.)

Maximum User Weight: 20 stone (125 Kg.)



Transfer Boards

Transfer Boards enable a patient to transfer and slide themselves from one level surface to another – bed to chair, commode or wheelchair.

- The boards are made from strong lacquered birch plywood with a wipe clean finish.

Long: 30" (76 cm.) with two handholds x Width 8 ¼" (21 cm)

Short: 24" (61 cm.) x Width 8 ¼" (21 cm)

Curved: 28 ¼ " (72 cm) x 10" (25 cm)

Maximum User Weight: 19 stone (120 Kg.)



Swivel Cushion

The Swivel Cushion is ideal in situations where help is needed with turning whilst sitting. Can be used on a chair, wheelchair, in bed and when a passenger in the car.

- Wipe clean plastic turntable which is fully rotational
- Washable fleece cover for comfort

Diameter: 16" (40.5 cm)

Maximum User Weight: 19 stone 9 lb. (125 Kg.)



Patient Handling UK,
2 East Court, Burley, Oakham,
Rutland LE15 0BD

Freephone: 0800 083 6071
Email: contact@patienthandling.org.uk
www.patienthandling.org.uk

Patient Handling UK

Suppliers of moving and handling aids for carers in the health industry and for individuals in their own homes.

Carefully designed products, constructed to high standards which meet the **CE Class 1 Medical Directives Standard.**

Orders receive prompt attention and dispatch is normally within 24 hours of order.

We are happy to deal with orders for a single product or for larger quantities.

Items can be bought via the website with secure online payment or orders can be emailed, faxed or posted.

Please contact us by Freephone or Email for advice on the suitability of a product or any other information required.



Sanitized® Slide Sheets

- Slippery silicon coated material with long lasting Sanitized® Antimicrobial protection.
- Very generous sizing for ease of use, folded or flat.
- Strong, resilient, low stretch and quiet to use.
- Write-on label for ward or patient marking.
- Washable. CE Class I Medical Device.

Standard Size - Red: 39" x 59" (1 m x 1.5 m)

Large Size - Aqua: 79" x 59" (2m x 1.5 m)



Maxi Tubular Slide Sheet

- Heavy duty roller slide sheet with strong webbing handles to both side.
- Ideal for the easier movement of heavier patients.
- Suitable for cross-bed and up-bed transfers.
- Can remain in place under a draw sheet or bed pad.
- Write-on label for patient or ward identification.
- Washable. CE Class I Medical Device.

Size: 28" x 61" (70 cm. x 153 cm.)



Walk & Transfer Belt

Assists when walking or during a transfer procedure.

- Padded belt with a non-slip lining to increase patient security.
- Easily adjustable buckle system to fit to waist.
- Four padded handles for carer assistance.
- Wipe clean finish. CE Class I Medical Device.

Small-red handles : Waist: 25" to 45" (63 cm.- 114 cm.)

Large-blue handles: Waist: 40" to 62" (100 cm.- 157 cm.)



One Way Seating Aid

Non-slip uni glide chair, bed or wheelchair cushion.

Prevents the user sliding forward or out of the seat.

Special internal fabric locks to minimize forward movement whilst allowing the user to be moved to the back of the chair.

Reduces the need for lap straps.

Washable.

CE Class I Medical Device.

One Size: 17" x 18" (42 cm. x 45 cm.)



Bed Ladder Strap

- Strong webbing strap with five handholds down the length.
- Assists user to sit upright in bed and to move from bed to chair.
- Provides exercise for the upper body during use.
- Fits to the bottom bed rail or the leg of most types of bed.
- Washable. CE Class I Medical Device.

Length: 75" (1.9 m.)



Stocking Slider

Handy nylon aid to fitting open toe compression stockings.

The Slider is fitted over the foot and heel.

The stocking is easily fitted over the Slider and up the leg.

The Slider is then removed by pulling on the loop handle.

One Size. Washable.



Freephone: 0800 083 6071

www.patienthandling.org.uk